


Disposable Blades

Low Profile

Description	Special Use
Shandon MX35 Premier + Low Profile Disposable Blade 80 x 8 x 0.25 mm 30° Cutting Angle	Sharper edge with increased durability. Provides improved performance for sectioning hard tissues and high-quality thin sections. Tests show 10-20% increase in blade life.
Shandon MX35 Premier Low Profile Disposable Blade 80 x 8 x 0.25 mm 35° Cutting Angle	Routine applications with a high proportion of harder tissue. Increased durability.
Shandon MB35 Premier Low Profile Disposable Blade 80 x 8 x 0.25 mm 35° Cutting Angle	Routine applications. Has consistent edge sharpness for high quality thin sections. Superior edge stability.
Shandon MB22 Premier Low Profile Disposable Blade 80 x 8 x 0.25 mm 22° Cutting Angle	Use for applications where extremely thin, high-quality sections are required such as renal biopsies and lymph nodes.




High Profile

Description	Special Use
Shandon HP35 High Profile Disposable Blade with PTFE Coating 75 x 14 x 0.31 mm 35° Cutting Angle	Routine applications with a high proportion of harder tissue. Increased durability with easier ribboning.
Shandon HP35n High Profile Disposable Blade without PTFE Coating 75 x 14 x 0.31 mm 35° Cutting Angle	Routine applications with a high proportion of harder tissue. Increased durability.

Structure of a Blade Edge


Superior Edge Profile

The shape of the blade edge is optimized to provide excellent sectioning performance, ensuring easier ribboning and improved sharpness.



Coating

Coating is optimized to provide durability and easy ribboning.



Disposable Blades Description	Packaging	Reorder No.
Shandon MX-35 Premier+ Low Profile 30° Cutting Angle	50/Pack	NEW 3052835
Shandon MX-35 Premier Low Profile 35° Cutting Angle	50/Pack	3051835
Shandon MB-35 Premier Low Profile 35° Cutting Angle	50/Pack	3050835
Shandon MB-22 Premier Low-Profile 22° Cutting Angle	50/Pack	3050822
Shandon HP-35 High Profile 35° Cutting Angle	50/Pack	3150734
Shandon HP-35n High Profile 35° Cutting Angle	50/Pack	3152735

