



TUBINGS - smooth / clear

reusable, autoclavable up to 134°C (273°F)

Material: Silicone (SI)

Box 1 roll = 25m (82ft.)

I.D. (mm)	O.D. (mm)	I.D. (inch)	O.D. (inch)	REF
1.0	2	.04"	.08"	68-010-05
1.0	3	.04"	.12"	68-010-10
1.5	3	.06"	.12"	68-015-075
2.0	4	.08"	.16"	68-020-10
3.0	5	.12"	.20"	68-030-10
3.0	6	.12"	.24"	68-030-15
4.0	6	.16"	.24"	68-040-10
4.0	7	.16"	.28"	68-040-15
4.0	8	.16"	.31"	68-040-20
5.0	7	.20"	.28"	68-050-10
5.0	8	.20"	.31"	68-050-15
5.0	9	.20"	.35"	68-050-20
6.0	8	.24"	.31"	68-060-10
6.0	9	.24"	.35"	68-060-15
6.0	10	.24"	.39"	68-060-20
6.0	12	.24"	.47"	68-060-30
7.0	10	.28"	.39	68-070-15
7.0	11	.28"	.43	68-070-20
7.0	13	.28"	.51"	68-070-30
8.0	12	.31"	.47"	68-080-20
8.0	14	.31"	.55"	68-080-30
9.0	13	.35"	.51"	68-090-20
10.0	14	.39"	.55"	68-100-20
10.0	16	.39"	.63"	68-100-30
12.0	18	.47"	.71"	68-120-30
12.0	20	.47"	.79"	68-120-40
14.0	20	.55"	.79"	68-140-30
15.0	21	.59"	.83"	68-150-30
16.0	22	.63"	.87"	68-160-30
18.0	24	.71"	.94"	68-180-30
20.0	24	.79"	.94"	68-200-20
20.0	26	.79"	1.02"	68-200-30
25.0	33	.98"	1.30"	68-250-40



Coil Extension Tubing without tube connectors

Material: Polyurethane (PUR), reusable, autoclavable up to 134°C (273°F)

Box 2.5m (8 1/4ft.)

Color	REF
red	20-10-442
blue	20-10-444
black	20-10-448

Box 1m (40"), max. 12.5m (41ft.) stretch length

Color	REF
red	20-20-442
blue	20-20-444
black	20-20-448

Dual Lumen Tubing, smooth, without tube connectors

Material: Silicone (SI), reusable, autoclavable up to 134°C (273°F)

Box 1m (40")

Color	REF
transparent 6mm O.D. / 3mm I.D. and 4mm O.D. / 2mm I.D.	68-360-24
transparent 8mm O.D. / 5mm I.D. and 4mm O.D. / 2mm I.D.	68-580-24

Material: Polyvinyl Chloride (PVC), single use

Box 1m (40")

Color	REF
red 6mm O.D. / 3mm I.D. and 6mm O.D. / 3mm I.D.	20-20-542
blue 6mm O.D. / 3mm I.D. and 6mm O.D. / 3mm I.D.	20-20-544

