

Pressure Infusion Cuffs

500ml and 1000ml

Manual system for rapid I.V. infusion of blood and solutions. Ideal for ambulances, operation rooms and intensive care stations.

500ml reusable

with hand inflator

REF

51-05-050

with hand inflator and manometer

51-06-050

with hand inflator and manometer
(shock resistant precision manometer)

51-08-050



REF 51-08-050

500ml

transparent, wrap around design, reusable REF

with hand inflator and pressure indicator

57-02-050

with hand inflator and manometer
(shock resistant precision manometer)

57-08-050



REF 57-08-050

1000ml

pocket design with zip, reusable

with hand inflator

REF

52-05-100

with hand inflator and manometer

52-06-100

with hand inflator and manometer
(shock resistant precision manometer)

52-08-100



REF 52-08-100

**1000ml
wrap around design, reusable**

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

53-05-100

53-06-100

53-08-100



**1000ml
for single patient use**

with female Luer Lock connection

with hand inflator and pressure indicator

with hand inflator

REF

56-01-100

56-02-100

56-03-100


**1000ml
transparent, wrap around design, reusable REF**

with hand inflator and pressure indicator

with hand inflator and manometer
(shock resistant precision manometer)

57-02-100

57-08-100



Pressure Infusion Cuffs

1500ml, 3000ml and 5000ml

For irrigation purposes e.g. during Arthroscopy, Laparoscopy and Hysteroscopy

1500ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer
(shock resistant precision manometer)

REF

52-01-150

52-08-150

52-30-150



3000ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer,
(shock resistant precision manometer)

REF

52-01-300

52-08-300

52-30-300

5000ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer,
(shock resistant precision manometer)

REF

52-01-500

52-08-500

52-30-500



Accessories for Pressure Infusion Cuffs



Manometer 300mmHg

shock resistant

REF 34-17-300



Manometer 300mmHg

REF 20-16-300



Manometer 300mmHg

with hand inflator and
release valve

REF 20-18-300



Manometer 300mmHg

with hand inflator, release valve
and male Luer Lock

REF 51-09-666



Hand Inflator

with release valve

REF 51-10-666



Handpump

with release valve and
male Luer Lock

REF 51-30-000



Foot Pump

with 2m (79") silicone hose,
release valve and male Luer Lock

REF 51-20-000



Hand Inflator

with pressure indicator, three-way
stop cock and male Luer Lock

REF 51-02-666



Hand Inflator

with three-way stop cock and
male Luer Lock

REF 51-03-666



Three-way stop cock

with two male and
one female Luer Lock

REF 51-02-000

